

## DNA Extraction

Before DNA extraction, phages should be at a sufficiently high concentration and should be filtered and free from bacteria.

### From Phages

- Add 10  $\mu$ L DNase I and 100  $\mu$ L of 10 x DNase buffer to up to 500  $\mu$ L of prepared phage sample and incubate at 37 °C for 1h
- Inactivate DNase I at 65 °C for 30 mins
- Add 40  $\mu$ L of proteinase K (20mg/mL) and incubate for 1h at 55 °C
- Incubate at -80 °C for 30 min
- Incubate at 70 °C for 30 min
- add up to 600uL of cold 100% ethanol to the tube and incubate overnight at -80 °C or -20°C
  - 100% ethanol should be kept at -20 °C before use
  - Sample can also be incubated at -20 °C
- Centrifuge the tubes at 16000-18000xg at 4 °C for 1h
- Remove the supernatant and resuspend the pellet in 100  $\mu$ L of TE buffer
- Purify the extracted DNA using ZYMO kit
  - Genomic DNA Clean & Concentrator
  - Use kit protocol as described
- Measure DNA concentration
  - If using nanodrop, record the ratios as well as the concentration
  - Qubit is more accurate