

# Phage Isolation

**Date:** 2025-11-06 ~ 2025-11-22 | **Operator:** Yuhang Gong | **Status:** Completed

## Objective

Isolate bacteriophages infecting EcAZ-1 (W strain) and EcAZ-2-OVA (R strain) from wastewater.

## Materials

### Bacterial Strains:

| Name       | Code | Source            | Notes          |
|------------|------|-------------------|----------------|
| EcAZ-1     | W    | Amir Zarrinar Lab | Wild type      |
| EcAZ-2-OVA | R    | Amir Zarrinar Lab | OVA-expressing |

**Reagents:** 2x LB, SM Buffer, 0.22um filter, Soft agar (0.7%), LB agar plates (1.5%)

**Equipment:** Centrifuge, 0.22um filtration apparatus, 37°C shaker, 37°C incubator, Micropipette

## Protocol

**Sample Processing:** Centrifuge wastewater → Filter supernatant (0.22um)

**Phage Enrichment:** Overnight host culture + 2xLB + wastewater supernatant → 37°C shaking overnight → Centrifuge & filter

**Screening & Purification:** Spot Test → Plaque Assay (double-layer) → Single clone isolation (5-6 rounds)

## Experimental Records

### Phase 1: Enrichment & Initial Screening (11-06 ~ 11-08)

| Date  | Operation  | Results                  |
|-------|--|--------------------------|
| 11-06 | Phage enrichment from wastewater; Spot test                            |                          |
| 11-07 | Spot test results; Plaque assay (10 <sup>-9</sup> ~10 <sup>-11</sup> ) | Phage presence confirmed |
| 11-08 | Plaque assay results; Single clone selection                           | R: 4 types; W: 2 types   |

### Phase 2: First Purification (11-09 ~ 11-10)

| Date  | Operation   | Results                                    |
|-------|---|--|
| 11-09 | Enrichment; Plaque assay (10 <sup>-9</sup> ~10 <sup>-12</sup> ) |  |
| 11-10 | Results   | R: failed (too concentrated); W: succeeded |

### Phase 3: Continued Purification (11-11 ~ 11-16)

| Date  | Operation  | Results                         |
|-------|--|---------------------------------|
| 11-11 | R: Plaque assay (10 <sup>-10</sup> , 10 <sup>-15</sup> ); W: Clone selection | R: still contaminated           |
| 11-12 | W: Clone selection; R: Filtered & stored                                     | W: Small plaques disappeared    |
| 11-13 | W: 4th purification; R: 3rd purification                                     |                                 |
| 11-14 | Post-purification enrichment   | W: 2 types; R: 3 types retained |
| 11-15 | W: 5th purification; R: 4th purification                                     |                                 |

| Date  | Operation           | Results  |
|-------|---------------------|--|
| 11-16 | Results; Enrichment | Essentially pure; Pipette tips > inoculation loops |

Phase 4: Final Purification (11-17 ~ 11-22)

| Date     | Operation                         | Results   |
|----------|-----------------------------------|---|
| 11-17~18 | Centrifuged; Prepared plates      |   |
| 11-19~20 | Plaque assay ( $10^{-10/12/14}$ ) | R3: need higher dilution; W1 split into W1 & W2 |
| 11-21~22 | Final plaque assay                | Pure stocks prepared                            |

Results

Successfully isolated 5 phage strains:

| Phage | Host       | Purification | Storage | Notes         |
|-------|------------|--------------|---------|---------------|
| R1    | EcAZ-2-OVA | 4-5 rounds   | 4°C     |               |
| R2    | EcAZ-2-OVA | 4-5 rounds   | 4°C     |               |
| R3    | EcAZ-2-OVA | 4-5 rounds   | 4°C     | High titer    |
| W1    | EcAZ-1     | 5-6 rounds   | 4°C     |               |
| W2    | EcAZ-1     | 5-6 rounds   | 4°C     | Split from W1 |

Conclusions

1. Isolated 5 phages from wastewater (R: 3 strains; W: 2 strains)
2. R3 has highest titer (requires  $>10^{-14}$  dilution)
3. W2 separated from W1 during purification
4. One unstable small-plaque phage was lost

Notes & Improvements

- High-titer phages require higher dilutions ( $10^{-14/16/18}$ )
- Use **pipette tips instead of inoculation loops** for clone picking
- Replace bacterial cultures weekly

Recorded: 2025-01-14